LEUKOCARE

formulation expertise meets data science

Analytical methods to support formulation development of Adeno-associated virus

Carina Rodenstein¹, Agatha Rosenthal², Bernd Schaefer², Natalia Markova², Andreas Seidl¹

¹: Leukocare AG, Munich, Germany, ²: Malvern Panalytical GmbH, Germany

Introduction

Background:


- Recombinant Adeno-associated viruses (rAAV) are widely used in gene therapy
- The benefit of optimal formulation is well recognized also for such drug products
- Analytical methods for characterization of rAAV are plentiful, typically cover capsid count, percentage of full rAAV particles, particle size, potential aggregate formation, capsid stability and infectivity.

Determination of relevant state-of-the-art analytical techniques for formulation development needs with AAV:

- Stability indicating & predictive power
- Low virus material requirement

Methods

Analytical methods typically used in AAV analysis and tested in this context:



High throughput option


Approach: Accelerated aging at 40 °C in reference buffer PBS/0.001% Poloxamer 188 (Ref) and in up to seven designed formulations (F01-F07). Model: AAV-2 and AAV-5 with insert Factor IX, a common genome size

Results

Particle integrity and aggregation: Analysis by AEX-HPLC, SE-HPLC, MALS, MADLS, and AUC as a standard


Analysis of AAV in reference buffer and formulations designed by Leukocare before & after thermal stress application (4 weeks at 40 °C):

- SE-HPLC & AUC indicated impact of thermal stress (reduction of full capsids) (AUC data not shown)
 - -> pattern of retained capsids after 4 weeks at 40 °C (stabilizing effect of the formulations) matched the pattern obtained with DSF (T_{m2}), nanoDSF, and DSC at t=0w
- SEC-MALS showed similar genome load as MADLS and qPCR-ELISA and additionally allowed aggregation analysis
- In AEX-HPLC, peaks of full & empty capsid were not separable after 4 weeks at 40°C, probably because of affected and changed net charge of capsids



A Full capsids (AAV-2) quantified by SE-HPLC/UV/MALS

Size distribution analysis of AAV-2 (SE-HPLC/UV) B



C Determination of genome load (AAV-5)

	Analytical method	Determined genome load [% full]		Expected genome load
		Mean	SD	[% full]
	MADLS	2.3**	n.a.	3.3
	SEC-MALS	3.9	0.1	(*qPCR+ELISA)
	AUC	6.0	n.a.	

* measurement performed by AAV vendor ** calculated based on qPCR results

SEC-MALS & AEX-HPLC successfully assessed full/empty capsids.

SE-HPLC /UV /MALS was stability indicating at accelerated aging condition, thereby discriminating between different formulations, and might replace AUC in formulation screening. SE-HPLC / UV / MALS combines genome load and aggregation analysis in one method.

Capsid protein unfolding and genome release: Analysis by nanoDSF, DSF, and DSC as a standard

Analysis of intrinsic tryptophane fluorescence by nanoDSF, before and after thermal stress application (4 weeks at 40 °C):

- Stability pattern (T_{on}, T_m) matched expectations based on experience with formulations
- No relevant impact of thermal stress observed by nanoDSF & DSC performance
- nanoDSF matched DSC results very well (tested on AAV-5)

Analysis by DSF / extrinsic dye-mediated DNA fluorescence (DSF / SYBR Gold)


- T_{m2} pattern matched expectations & nanoDSF results on protein unfolding
- T_{m1} indicated genome release prior to capsid protein unfolding; despite high capsid protein stability, the pattern differed from T_{m2} results



All methods (NanoDSF, DSF, DSC) gave comparable results, could discriminate formulations, & predicted stability (t=0w vs SEC-MALS t=4w_40°C). But these methods did not show relevant effects of thermal stress (t=0w vs t=4w_40°C, not stability indicating), which were detected by SEC-MALS.

DSF additionally provided information on DNA release prior protein unfolding detection, however, requires a fluorescent dye.

Summary



contact@leukocare.com

Conclusion

For each relevant formulation parameter at least one suitable method could be identified that fulfilled our requirements. These were defined to be stability indicating or predictive.

Acknowledgements

• NTA

Excellent support at Leukocare: Technical work in laboratory

Viviane Firmery, Christiane Langenegger Joshua Schäfer, Michael Sporer

Scientific consultation:

Eva Reinauer



We gratefully thank Malvern Panalytical for providing the following equipment MALS Detector DSC

MADLS

We gratefully thank Protagen, Germany, for performing the AUC measurements